



Welcome to MedConnect's Document Viewer / Editor

To begin, click the  icon to upload an existing document located on your computer.

If you have a scanner installed, click the  icon to select your scanner and then click the  icon to upload an image directly from the scanner.

After your image has been loaded, edit the document using any of the tools provided. Refer to the list below for a description of the available tools.

Tools available:

File Operations:

	Opens a file from your computer
	Scans in a new batch of images
	Prepares the images for printing from your browser
	Downloads the current file to your computer

Navigation:

	Moves to the previous image
	Shows the current image number; Typing a number navigates to that image
	Moves to the next image

Paging:

	Deletes the current image
	Inserts a new image
	Scans a new image and replaces the current image

Zoom:

	Zoom In
	Zoom Out
	Fit the image to the width of the viewer
	Fit the entire image to the viewer
	View the image at its real size

Processing:

	Rotates the image clockwise 90 degrees
	Rotates the image counter-clockwise 90 degrees
	Rotates the image 180 degrees
	Crops the image using a selection rectangle
	Removes specks from the image
	Deskews the images
	Reduces the color image to black and white
	Removes hole punches from the image
	Remove the black border from the image