

Patient Information	Specimen Information	Client Information
HIGDON, LEON DOB: 03/20/1966 AGE: 54 Gender: M Fasting: Y Phone: 251.622.1918 Patient ID: NG	Specimen: AL280018Q Requisition: 6028675 Collected: 05/21/2020 / 08:45 CDT Received: 05/22/2020 / 02:12 CDT Reported: 05/22/2020 / 08:38 CDT	Client #: 49100638 QATL0000 TUCK, PAMELA P AUM NURSING CARE CENTER 7461 EAST DRIVE #102 MONTGOMERY, AL 36117

Test Name	In Range	Out Of Range	Reference Range	Lab
LIPID PANEL, STANDARD				
CHOLESTEROL, TOTAL		208 H	<200 mg/dL	AT
HDL CHOLESTEROL	55		> OR = 40 mg/dL	AT
TRIGLYCERIDES		221 H	<150 mg/dL	AT

If a non-fasting specimen was collected, consider repeat triglyceride testing on a fasting specimen if clinically indicated.
Jacobson et al. J. of Clin. Lipidol. 2015;9:129-169.

LDL-CHOLESTEROL	120 H	mg/dL (calc)	AT
Reference range: <100			

Desirable range <100 mg/dL for primary prevention;
<70 mg/dL for patients with CHD or diabetic patients with > or = 2 CHD risk factors.

LDL-C is now calculated using the Martin-Hopkins calculation, which is a validated novel method providing better accuracy than the Friedewald equation in the estimation of LDL-C.

Martin SS et al. JAMA. 2013;310(19): 2061-2068
(<http://education.QuestDiagnostics.com/faq/FAQ164>)

CHOL/HDL-C RATIO	3.8	<5.0 (calc)	AT
NON HDL CHOLESTEROL	153 H	<130 mg/dL (calc)	AT

For patients with diabetes plus 1 major ASCVD risk factor, treating to a non-HDL-C goal of <100 mg/dL (LDL-C of <70 mg/dL) is considered a therapeutic option.

COMPREHENSIVE METABOLIC PANEL			AT
GLUCOSE	101 H	65-99 mg/dL	

Fasting reference interval

For someone without known diabetes, a glucose value between 100 and 125 mg/dL is consistent with prediabetes and should be confirmed with a follow-up test.

UREA NITROGEN (BUN)	12	7-25 mg/dL
CREATININE	1.08	0.70-1.33 mg/dL

For patients >49 years of age, the reference limit for Creatinine is approximately 13% higher for people identified as African-American.

eGFR NON-AFR. AMERICAN	77	> OR = 60 mL/min/1.73m2
eGFR AFRICAN AMERICAN	90	> OR = 60 mL/min/1.73m2
BUN/CREATININE RATIO	NOT APPLICABLE	6-22 (calc)
SODIUM	138	135-146 mmol/L
POTASSIUM	4.3	3.5-5.3 mmol/L
CHLORIDE	102	98-110 mmol/L
CARBON DIOXIDE	25	20-32 mmol/L

Patient Information	Specimen Information	Client Information
HIGDON, LEON DOB: 03/20/1966 AGE: 54 Gender: M Fasting: Y Patient ID: NG	Specimen: AL280018Q Collected: 05/21/2020 / 08:45 CDT Received: 05/22/2020 / 02:12 CDT Reported: 05/22/2020 / 08:38 CDT	Client #: 49100638 TUCK, PAMELA P

Test Name	In Range	Out Of Range	Reference Range	Lab
CALCIUM	9.9		8.6-10.3 mg/dL	
PROTEIN, TOTAL	7.3		6.1-8.1 g/dL	
ALBUMIN	4.6		3.6-5.1 g/dL	
GLOBULIN	2.7		1.9-3.7 g/dL (calc)	
ALBUMIN/GLOBULIN RATIO	1.7		1.0-2.5 (calc)	
BILIRUBIN, TOTAL	0.6		0.2-1.2 mg/dL	
ALKALINE PHOSPHATASE	69		35-144 U/L	
AST	20		10-35 U/L	
ALT	27		9-46 U/L	
TSH W/REFLEX TO FT4	1.84		0.40-4.50 mIU/L	AT
CBC (INCLUDES DIFF/PLT)				AT
WHITE BLOOD CELL COUNT	5.8		3.8-10.8 Thousand/uL	
RED BLOOD CELL COUNT	5.14		4.20-5.80 Million/uL	
HEMOGLOBIN	16.4		13.2-17.1 g/dL	
HEMATOCRIT	47.2		38.5-50.0 %	
MCV	91.8		80.0-100.0 fL	
MCH	31.9		27.0-33.0 pg	
MCHC	34.7		32.0-36.0 g/dL	
RDW	12.5		11.0-15.0 %	
PLATELET COUNT	352		140-400 Thousand/uL	
MPV	9.5		7.5-12.5 fL	
ABSOLUTE NEUTROPHILS	3584		1500-7800 cells/uL	
ABSOLUTE LYMPHOCYTES	1357		850-3900 cells/uL	
ABSOLUTE MONOCYTES	545		200-950 cells/uL	
ABSOLUTE EOSINOPHILS	232		15-500 cells/uL	
ABSOLUTE BASOPHILS	81		0-200 cells/uL	
NEUTROPHILS	61.8		%	
LYMPHOCYTES	23.4		%	
MONOCYTES	9.4		%	
EOSINOPHILS	4.0		%	
BASOPHILS	1.4		%	

Patient Information	Specimen Information	Client Information
HIGDON, LEON DOB: 03/20/1966 AGE: 54 Gender: M Fasting: Y Patient ID: NG	Specimen: AL280018Q Collected: 05/21/2020 / 08:45 CDT Received: 05/22/2020 / 02:12 CDT Reported: 05/22/2020 / 08:38 CDT	Client #: 49100638 TUCK, PAMELA P

Endocrinology

Test Name	Result	Reference Range	Lab
VITAMIN D,25-OH,TOTAL,IA	26 L	30-100 ng/mL	AT
Vitamin D Status 25-OH Vitamin D: Deficiency: <20 ng/mL Insufficiency: 20 - 29 ng/mL Optimal: > or = 30 ng/mL For 25-OH Vitamin D testing on patients on D2-supplementation and patients for whom quantitation of D2 and D3 fractions is required, the QuestAssureD(TM) 25-OH VIT D, (D2,D3), LC/MS/MS is recommended: order code 92888 (patients >2yrs).			
Physician Comments:			

End Notes:

VITAMIN D,25-OH,TOTAL,IA

AT

For additional information, please refer to <http://education.QuestDiagnostics.com/faq/FAQ199> (This link is being provided for informational/ educational purposes only.)

PERFORMING SITE:

AT QUEST DIAGNOSTICS-ATLANTA, 1777 MONTREAL CIRCLE, TUCKER, GA 30084-6802 Laboratory Director: ANDREW N YOUNG,MD,PHD, CLIA: 11D0255931

LIST OF RESULTS PRINTED IN THE OUT OF RANGE COLUMN:

CHOLESTEROL, TOTAL	208 H	<200 mg/dL	AT
TRIGLYCERIDES	221 H	<150 mg/dL	AT

If a non-fasting specimen was collected, consider repeat triglyceride testing on a fasting specimen if clinically indicated.
Jacobson et al. J. of Clin. Lipidol. 2015;9:129-169.

LDL-CHOLESTEROL	120 H	mg/dL (calc)	AT
Reference range: <100			

Desirable range <100 mg/dL for primary prevention;
<70 mg/dL for patients with CHD or diabetic patients with > or = 2 CHD risk factors.

LDL-C is now calculated using the Martin-Hopkins calculation, which is a validated novel method providing better accuracy than the Friedewald equation in the estimation of LDL-C.

Martin SS et al. JAMA. 2013;310(19): 2061-2068
(<http://education.QuestDiagnostics.com/faq/FAQ164>)

NON HDL CHOLESTEROL	153 H	<130 mg/dL (calc)	AT
---------------------	-------	-------------------	----

For patients with diabetes plus 1 major ASCVD risk factor, treating to a non-HDL-C goal of <100 mg/dL (LDL-C of <70 mg/dL) is considered a therapeutic option.

GLUCOSE	101 H	65-99 mg/dL	AT
---------	-------	-------------	----

Fasting reference interval

Patient Information	Specimen Information	Client Information
HIGDON, LEON DOB: 03/20/1966 AGE: 54 Gender: M Fasting: Y Patient ID: NG	Specimen: AL280018Q Collected: 05/21/2020 / 08:45 CDT Received: 05/22/2020 / 02:12 CDT Reported: 05/22/2020 / 08:38 CDT	Client #: 49100638 TUCK, PAMELA P

LIST OF RESULTS PRINTED IN THE OUT OF RANGE COLUMN:

For someone without known diabetes, a glucose value between 100 and 125 mg/dL is consistent with prediabetes and should be confirmed with a follow-up test.

VITAMIN D, 25-OH, TOTAL, IA	26 L	30-100 ng/mL	AT
Vitamin D Status	25-OH Vitamin D:		
Deficiency:	<20 ng/mL		
Insufficiency:	20 - 29 ng/mL		
Optimal:	> or = 30 ng/mL		

For 25-OH Vitamin D testing on patients on D2-supplementation and patients for whom quantitation of D2 and D3 fractions is required, the QuestAssured(TM) 25-OH VIT D, (D2,D3), LC/MS/MS is recommended: order code 92888 (patients >2yrs).